




New gel size




New sample size




# Sample view and fragment view

Click a microtube tab or its lane to select the sample and bring up the sample view.  
Right-click on lane to make a band.  
Command-click to create a loupe  
Escape/Delete to close a selected loupe..

- TODO:
- ✓ New button
  - ✓ Global settings
    - ✓ Toggle color bands
    - ✓ Toggle mouse => bp rollover
  - ✓ Gel settings
  - ✓ Notate all click key combos.
  - Smooth and discrete multimer transitions
  - Monomer default is clearer.
  - (Talk to Shawn)
  - ✓ Electrical lines
  - ✓ N lanes? — just show a demo with shrunken bits.




Running buffer +  
electrical lines




# More lanes

If we add a slider for # of lanes, then we just  
resize the lanes as needed.




## Sample buffer and dyes

**Click background** to bring up dyes; **right-click** brings up buffer + dyes (Or we just simplify and show all settings together on background click).



**Command-Plus** or **Command-N** makes a new fragment (or button). **Delete/Backspace** to delete a fragment.

# Gel settings



Push app settings icon to toggle app settings menu.

Push gel settings icon to toggle gel settings: time, voltage, well damage, and rotation.

# Gel and buffer settings

Loupe

Gel

Buffer

.5x TBE

1x TAE

Tris

Borate


Acetate

EDTA

Right click gel settings to get buffer parameters, too.  
(Or we just show them all together all the time.)

# Icon exploration

## Settings



37x37

32x32

## New species



32x32



Loupe






New button  
Global settings  
Gel settings



Notate all click key combos.




 Sample




Species



Loupe



 Sample







Buffer + Dyes

Buffer    .5x TBE    1x TAE

Tris	<div></div>	
Borate	<div></div>	
Acetate	<div></div>	
EDTA	<div></div>	

Dyes

Xylene cyanol	<div></div>	
Cresol red	<div></div>	
Bromphenol blue	<div></div>	
Orange G	<div></div>	

Loupe

